


# Axle<sup>®</sup>




Interspinous Fusion System



Modular Interspinous Fusion Device


# Axle<sup>®</sup>

## Fusion System



The Axle interspinous fusion system offers multiple configurations allowing the surgeon to precisely tailor the implant to the patient's anatomy. Implanted via a tissue sparing approach, Axle presents a viable option for surgeons looking for a minimally invasive alternative to pedicle screws for rigid posterior fixation.

- Clinically proven alternative to bilateral pedicle screw fixation in a retrospective TLIF case series while potentially offering reduced operating room time and hospital stay.\*
- Modular systems enable a precise fit to patient anatomy with either PEEK or titanium inserts to aid in spinous process distraction. Inserts can be loaded with bone graft to enhance fusion potential.
- Comprehensive system including flared plate option to capture L5-S1 as well as complementary allograft insert to enhance spinous process distraction.



### Also available:



### H Graft<sup>®</sup>

Interspinous Allograft

H-shaped cortical allograft block available in 2mm increments from 8mm to 20mm.



☎ 888.886.9354  
✉ [cs@xtantmedical.com](mailto:cs@xtantmedical.com)  
🌐 [xtantmedical.com](http://xtantmedical.com)  
📍 664 Cruiser Ln,  
Belgrade, MT 59714

INDICATIONS: See Package Insert for a more complete listing of indications, contraindications, warnings, precautions, and other important information.  
LIMITED WARRANTY and DISCLAIMER: X-spine<sup>®</sup> Systems, Inc. products have a limited warranty against defects and workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are disclaimed. WARNING: In the USA, this product has labeling limitations. See package insert for complete information. CAUTION: USA Law restricts these devices to sale by or on the order of a physician.  
\*Data on file at Xtant Medical.

Axle is a registered trademark of X-Spine Systems, a subsidiary of Xtant Medical.  
© Xtant Medical. All Rights Reserved. FM-C-MRK-27 (E) 06/2020